

# Schedule of Accreditation

issued by

## United Kingdom Accreditation Service

21 - 47 High Street, Feltham, Middlesex, TW13 4UN, UK

 <p><b>UKAS</b> CALIBRATION</p> <p>0698</p> <p>Accredited to ISO/IEC 17025:2005</p>	<p><b>Biohit Ltd trading as Pipette Doctor</b></p> <p>Issue No: 008    Issue date: 23 July 2010</p>	
	<p>Unit 1 Barton Hill Way Torquay Devon TQ2 8JG</p>	<p>Contact: Mr Bernie Atkinson Tel: +44 (0)1803 315900 Fax: +44 (0)1803 313584 E-Mail: <a href="mailto:berniea@pipettedoctor.com">berniea@pipettedoctor.com</a> Website: <a href="http://www.pipettedoctor.com">www.pipettedoctor.com</a></p>
<p>Calibration performed by the Organisations at the locations specified below</p>		

### Locations covered by the organisation and their relevant activities

#### Laboratory locations:

Location details	Activity	Location code
<p><b>Address</b> Unit 1 Barton Hill Way Torquay Devon TQ2 8JG</p> <p><b>Local contact</b> Mr B Atkinson</p>	Volume	A

#### Site activities performed away from the locations listed above:

Location details	Activity	Location code
<p>Any premises The customers' site or premises must be suitable for the nature of the particular calibrations undertaken and will be the subject of contract review arrangements between the laboratory and the customer.</p> <p><b>Contact:</b> Mr B Atkinson</p>	Volume	B



0698  
Accredited to  
ISO/IEC 17025:2005

**Schedule of Accreditation**  
issued by  
**United Kingdom Accreditation Service**  
21 - 47 High Street, Feltham, Middlesex, TW13 4UN, UK

**Biohit Ltd trading as Pipette Doctor**  
**Issue No: 008 Issue date: 23 July 2010**

Calibration performed by the Organisation at the locations specified

DETAIL OF ACCREDITATION

Measured Quantity Instrument or Gauge	Range	Calibration and Measurement Capability Expressed as an Expanded Uncertainty ( $k=2$ )	Remarks	Location Code
VOLUME of liquids (See Note 1)	Nominal value		Note 1. For water delivered from piston and/or plunger operated volumetric apparatus, using procedures agreed with UKAS.  Note 2. The CMCs quoted apply to calibrations performed with either 10 aliquots or 5 aliquots taken from each pipette at the time of the calibration.	A
	From 0.2 $\mu$ L to 10.0 $\mu$ L Above 10.0 $\mu$ L to 20.0 $\mu$ L Above 20.0 $\mu$ L to 100 $\mu$ L Above 100 $\mu$ L to 200 $\mu$ L Above 200 $\mu$ L to 500 $\mu$ L Above 0.5 mL to 1 mL Above 1 mL to 2 mL Above 2 mL to 5 mL Above 5 mL to 10 mL	0.15 $\mu$ L 0.25 $\mu$ L 0.40 $\mu$ L 0.60 $\mu$ L 1.4 $\mu$ L 0.003 5 mL 0.008 0 mL 0.020 mL 0.030 mL		
	Above 10.0 $\mu$ L to 20.0 $\mu$ L Above 20.0 $\mu$ L to 100 $\mu$ L Above 100 $\mu$ L to 200 $\mu$ L Above 200 $\mu$ L to 500 $\mu$ L Above 0.5 mL to 1 mL Above 1 mL to 2 mL	0.40 $\mu$ L 0.60 $\mu$ L 0.80 $\mu$ L 1.5 $\mu$ L 0.004 0 mL 0.007 0 mL		B
END				